

This document is prepared by: NHS National Procurement Enquiries to: **nss.pcfparamed@nhs.scot** The following amendments to Stoma Appliances take effect from 1st December 2024.

Please note that potential additions to the Scottish Stoma Appliances Reimbursement List must be referred through National Procurement (NSS)

ADDITIONS

▲ These products qualify for the stoma customisation fee

COLOSTOMY BAGS

B. Braun Medical Ltd	Softima Active Convex Closed 1 piece closed pouch with soft convex skin protector, filter and protective split cover Midi Beige split cover					
	Cut-to-fit	15-45mm	4612716E	10	3450	
SKIN FILLERS AND PROTECTIVES						
Teva UK Ltd	Sudocrem Care & Protect AssuraDerm Ointment			100g	450	

TWO PIECE OSTOMY SYSTEMS

Coloplast Ltd

▲ SenSura Mio Flex 2-piece baseplate with integral soft convexity, BodyFit Technology and adhesive coupling					
35mm coupling, starter hole 10-20mm	16880	5	2190		
50mm coupling, starter hole 15-30mm	16882	5	2190		
50mm coupling, starter hole 15-40mm	16886	5	2190		
70mm coupling, starter hole 15-53mm	16889	5	2190		

DELETIONS

TWO PIECE OSTOMY SYSTEMS

Peak Medical Ltd

Combimate Anatomical Closed pouch, Beige non-woven both sides, high capacity filter

Medium/Large	57mm flange	CMLB257	30	4750
Combimate SuperFlex 2	2-piece Flange			
	45mm flange			
	pre-cut 22mm	FFB4522	5	2138
	pre-cut 25mm	FFB4525	5	2138
	pre-cut 28mm	FFB4528	5	2138
	pre-cut 32mm	FFB4532	5	2138
	pre-cut 35mm	FFB4535	5	2138
57mm flange				
	pre-cut 32mm	FFB5732	5	2138
	pre-cut 35mm	FFB5735	5	2138
	pre-cut 38mm	FFB5738	5	2138

Combimate Supersoft flange						
45mm flange						
pre-cut 22mm	FSS4522	5	2138			
pre-cut 25mm	FSS4525	5	2138			
pre-cut 28mm	FSS4528	5	2138			
pre-cut 32mm	FSS4532	5	2138			
pre-cut 35mm	FSS4535	5	2138			
57mm flange						
13mm starter hole	FSS5713	5	2138			
pre-cut 32mm	FSS5732	5	2138			
pre-cut 35mm	FSS5735	5	2138			
pre-cut 38mm	FSS5738	5	2138			
70mm flange						
13mm starter hole	FSS7013	5	2138			
Combimate Supersoft flange with shallow convexity						
45mm flange						
pre-cut 22mm	FC4522	5	2458			
pre-cut 25mm	FC4525	5	2458			
pre-cut 28mm	FC4528	5	2458			
57mm flange						
pre-cut 32mm	FC5732	5	2458			
pre-cut 35mm	FC5735	5	2458			
pre-cut 38mm	FC5738	5	2458			
pre-cut 41mm	FC5741	5	2458			
Combimate Supersoft flange with deep convexity						
45mm flange						
pre-cut 19mm	FXC4519	5	2458			
pre-cut 22mm	FXC4522	5	2458			
57mm flange						
pre-cut 25mm	FXC5725	5	2458			
pre-cut 28mm	FXC5728	5	2458			
pre-cut 32mm	FXC5732	5	2458			

AMENDMENTS

 \star These product codes will change in the December 2024 Drug Tariff

These products must be dispensed with a supply of wipes and disposal bags

▲ These products qualify for the stoma customisation fee

These product codes will change in the March 2025 Drug Tariff

IRRIGATION/WASH - OUT APPLIANCES

CliniSupplies Ltd (formerly Renew Medical UK)

Now Reads

CliniSupplies Ltd

TWO PIECE OSTOMY SYSTEMS

Peak Medical Ltd

★ ■ Combimate Closed pouch, Beige non-woven both sides, high capacity filter

Mediu	IM			
	45mm flange	CMB245	30	4750
Large				
	45mm flange	CLB245	30	4750
	57mm flange	CLB257	30	4750
	70mm flange	CLB270	30	4750

	ch, Beige non-woven both sides, hi 1edium	gh capacity filter		
	45mm flange	421995	30	4750
La	arge			
	45mm flange	421991	30	4750
	57mm flange	421992	30	4750
	70mm flange	421989	30	4750
	Closed pouch, Beige non-woven be ledium/Large	oth sides, high capacity filter		
*	45mm flange	CMLB245	30	4750
Now Reads				
Μ	1edium/Large 45mm flange	421998	30	4750
▲ Combimate oval soft co	nnvex flange			
	mm flange			
~ J7	starter hole 13mm	FSC05713-44	5	2408
Now Reads		13003/13 44	5	2400
	7mm flange			
•	starter hole 13mm	422785	5	2408
			-	
★ ▲ Combimate soft conv	vex flange			
57	7mm flange			
	starter hole 13mm	FSC5713-36	5	2438
Now Reads				
Combimate soft convex	-			
57	7mm flange		_	
	starter hole 13mm	422777	5	2438
Transparent	bouch, with non-return valve and so	oft bottom outlet tap		
Medium			20	6405
	5mm flange	UMT245B	20	6105
	7mm flange	UMT257B	20	6105
Now Reads	mm flange	UMT245B	20	6105
	'mm flange	UMT257B	20	6105
* 57	IIIII lange	010112378	20	0103
▲ Combimate SuperFlex 2	2-piece Flange			
•	5mm flange			
*	starter hole 13mm	FFB4513	5	2138
57	7mm flange			
*	starter hole 13mm	FFB5713	5	2138
Now Reads				
45	5mm flange			
	starter hole 13mm	422020	5	2138
57	7mm flange			
	starter hole 13mm	422024	5	2138
▲ Combimate Supersoft f	lange			
•	5mm flange			
	13mm starter hole	FSS4513	5	2138
		r334513	5	2130
Now Reads	Creare flammer			
45	5mm flange	400007	-	2420
	13mm starter hole	422037	5	2138

	0. On a second state of the second state of th					
Combinate Supersoft flange with shallow convexity						
	45mm flange		504540	-	2450	
	\star 13mm starter hole		FC4513	5	2458	
	57mm flange		505740	_	2450	
	\star 13mm starter hole		FC5713	5	2458	
Now Reads						
	45mm flange			_		
	13mm starter hole		422029	5	2458	
	57mm flange					
	13mm starter hole		422031	5	2458	
🔺 Combimate Superso	oft flange with deep conv	vexity				
	45mm flange					
	\star 13mm starter hole		FXC4513	5	2458	
	57mm flange					
	★ 13mm starter hole		FXC5713	5	2458	
Now Reads						
	45mm flange					
	13mm starter hole		422042	5	2458	
	57mm flange					
	13mm starter hole		422044	5	2458	
UROSTOMY BAGS						
Peak Medical Ltd						
Uromate 1 piece extra outlet tap	small urostomy pouch, s	emi-transparent, beige n	on-woven, non-returr	n valve, so	oft bottom	
		10mm starter hole	UXST110B	20	11757	
Uromate 1 piece ext	ra small urostomy pouch	n, semi-transparent, beig				
outlet tap		,			,	
		10mm starter hole	UXST110B	20	11757	
Uromate 1-niece urost	omy pouch with split cov	ver, non-return valve, 4-p	oint helt attachment a	and soft k	ottom outlet	
tap		ei, non-return valve, 4-μ			ottom outlet	
	Beige					
	Large					
	starter hole 13-35mm		ULTP113B-35	20	11612	
	13-44mm		ULTP113B-44	20	11612	
	13-54mm		ULTP113B-54	20	11612	
Now Reads						
Uromate 1-piece urostomy pouch with split cover, non-return valve, 4-point belt attachment and soft bottom						
outlet tap						
	Beige					
	Large					
	starter hole 13-35mm		ULTP113B-35	20	11612	
	13-44mm		ULTP113B-44	20	11612	
	13-54mm		ULTP113B-54	20	11612	